Rekommenderade metoder

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Referenser

1. World Health Organization. Laboratory biosafety manual. Third edition. Geneva, WHO. 2004.
<https://www.who.int/csr/resources/publications/biosafety/Biosafety7.pdf?ua=1>
2. Principles and procedures for blood cultures; approved standard—1st ed. CLSI document M47-A.: Clinical and Laboratory Standards Institute, Wayne, PA., 2007.
3. Working Party of the British Society for Antimicrobial Chemotherapy. Guidelines for the diagnosis and antibiotic treatment of endocarditis in adults: a report of the Working Party of the British Society for Antimicrobial Chemotherapy. *J Antimicrob Chemother* 67: British Society for Antimicrobial Chemotherapy, 2012; 269-89.
4. Public Health England. (2017). Investigation of intravascular cannulae and associated specimens. UK Standards for Microbiology Investigations. B 20 Issue 6.1. London: National Health Service (NHS), 2017; 1-24.
<https://www.gov.uk/government/collections/standards-for-microbiology-investigations-smi>
5. Public Health England. (2019). Investigation of blood cultures (for organisms other than Mycobacterium species). UK Standards for Microbiology Investigations. B 37 Issue 8.2. London: National Health Service (NHS), 2019; 1-51.
<https://www.gov.uk/government/collections/standards-for-microbiology-investigations-smi>
6. Vårdprogram, Infektiös endokardit, rev 2021. Svenska infektionsläkarföreningen, 2021.
7. Acuña M, O’Ryan M, Cofré J et al. Differential time to positivity and quantitative cultures for noninvasive diagnosis of catheter-related blood stream infection in children. *The Pediatric infectious disease journal* 2008; **27**: 681-5.
8. Alestig K, Hogevik H, Olaison L. Infective endocarditis: a diagnostic and therapeutic challenge for the new millennium. *Scand J Infect Dis* 2000; **32**: 343-56.
9. Almuhayawi M, Altun O, Abdulmajeed AD et al. The Performance of the Four Anaerobic Blood Culture Bottles BacT/ALERT-FN, -FN Plus, BACTEC-Plus and -Lytic in Detection of Anaerobic Bacteria and Identification by Direct MALDI-TOF MS. *PloS one* 2015; **10**: e0142398.
10. Altun O, Almuhayawi M, Ullberg M et al. Clinical evaluation of the FilmArray blood culture identification panel in identification of bacteria and yeasts from positive blood culture bottles. *J Clin Microbiol* 2013; **51**: 4130-6.
11. Altun O, Botero-Kleiven S, Carlsson S et al. Rapid identification of bacteria from positive blood culture bottles by MALDI-TOF MS following short-term incubation on solid media. *J Med Microbiol* 2015; **64**: 1346-52.
12. Arbetsmiljöverket. Smittrisker (AFS 2018:4). <https://www.av.se>
13. Arendrup M, Jensen IP, Justesen T. Diagnosing bacteremia at a Danish hospital using one early large blood volume for culture. *Scand J Infect Dis* 1996; **28**: 609-14.
14. Baddour LM, Epstein AE, Erickson CC et al. Update on cardiovascular implantable electronic device infections and their management: a scientific statement from the American Heart Association. *Circulation* 2010; **121**: 458-77.
15. Baddour LM, Epstein AE, Erickson CC et al. A summary of the update on cardiovascular implantable electronic device infections and their management: a scientific statement from the American Heart Association. *J Am Dent Assoc* 2011; **142**: 159-65.
16. Banerjee R, Ozenci V, Patel R. Individualized Approaches Are Needed for Optimized Blood Cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2016; **63**: 1332-9.

17. Barenfanger J, Graham DR, Kolluri L et al. Decreased mortality associated with prompt Gram staining of blood cultures. *American journal of clinical pathology* 2008; **130**: 870-6.
18. Baron EJ, Miller JM, Weinstein MP et al. A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2013 recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM)(a). *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2013; **57**: e22-e121.
19. Baron EJ, Scott JD, Tompkins LS. Prolonged incubation and extensive subculturing do not increase recovery of clinically significant microorganisms from standard automated blood cultures. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2005; **41**: 1677-80.
20. Beekmann SE, Diekema DJ, Chapin KC et al. Effects of rapid detection of bloodstream infections on length of hospitalization and hospital charges. *J Clin Microbiol* 2003; **41**: 3119-25.
21. Ben-Ami R, Weinberger M, Orni-Wasserlauff R et al. Time to blood culture positivity as a marker for catheter-related candidemia. *J Clin Microbiol* 2008; **46**: 2222-6.
22. Binkhamis K, Forward K. Effect of the initial specimen diversion technique on blood culture contamination rates. *J Clin Microbiol* 2014; **52**: 980-1.
23. Bouza E, Alcalá L, Muñoz P et al. Can microbiologists help to assess catheter involvement in candidaemic patients before removal? *Clin Microbiol Infect* 2013; **19**: E129-35.
24. Bouza E, Alvarado N, Alcalá L et al. A prospective, randomized, and comparative study of 3 different methods for the diagnosis of intravascular catheter colonization. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2005; **40**: 1096-100.
25. Bouza E, Burillo A, Muñoz P. Catheter-related infections: diagnosis and intravascular treatment. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases* 2002; **8**: 265-74.
26. Bouza E, Burillo A, Muñoz P et al. Mixed bloodstream infections involving bacteria and *Candida* spp. *The Journal of antimicrobial chemotherapy* 2013; **68**: 1881-8.
27. Bouzidi H, Emirian A, Marty A et al. Differential time to positivity of central and peripheral blood cultures is inaccurate for the diagnosis of *Staphylococcus aureus* long-term catheter-related sepsis. *J Hosp Infect* 2018; **99**: 192-9.
28. Brun-Buisson C, Abrouk F, Legrand P et al. Diagnosis of central venous catheter-related sepsis. Critical level of quantitative tip cultures. *Arch Intern Med* 1987; **147**: 873-7.
29. Calabrese F, Carturan E, Thiene G. Cardiac infections: focus on molecular diagnosis. *Cardiovasc Pathol* 2010; **19**: 171-82.
30. Cardoso CL, Muraro CB, Siqueira VL et al. Simplified technique for detection of significant bacteriuria by microscopic examination of urine. *J Clin Microbiol* 1998; **36**: 820-3.
31. Cascone VJ, Cohen RS, Dodson NP et al. Implications of culture collection after the first antimicrobial dose in septic emergency department patients. *The American journal of emergency medicine* 2019; **37**: 947-51.
32. Catton JA, Dobbins BM, Kite P et al. In situ diagnosis of intravascular catheter-related bloodstream infection: a comparison of quantitative culture, differential time to positivity, and endoluminal brushing. *Critical care medicine* 2005; **33**: 787-91.
33. Chaiyakunapruk N, Veenstra DL, Lipsky BA et al. Vascular catheter site care: the clinical and economic benefits of chlorhexidine gluconate compared with povidone

- iodine. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2003; **37**: 764-71.
34. Chen IH, Nicolau DP, Kuti JL. Recovery of Gram-Negative Bacteria from Aerobic Blood Culture Bottles Containing Antibiotic Binding Resins after Exposure to beta-Lactam and Fluoroquinolone Concentrations. *J Clin Microbiol* 2019; **57**.
 35. Chen IH, Nicolau DP, Kuti JL. Effect of Clinically Meaningful Antibiotic Concentrations on Recovery of Escherichia coli and Klebsiella pneumoniae Isolates from Anaerobic Blood Culture Bottles with and without Antibiotic Binding Resins. *J Clin Microbiol* 2019; **57**.
 36. Cheng MP, Stenstrom R, Paquette K et al. Blood Culture Results Before and After Antimicrobial Administration in Patients With Severe Manifestations of Sepsis: A Diagnostic Study. *Ann Intern Med* 2019; **171**: 547-54.
 37. Chung Y, Kim IH, Han M et al. A comparative evaluation of BACT/ALERT FA PLUS and FN PLUS blood culture bottles and BD BACTEC Plus Aerobic and Anaerobic blood culture bottles for antimicrobial neutralization. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology* 2019; **38**: 2229-33.
 38. Clancy CJ, Nguyen MH. Finding the "missing 50%" of invasive candidiasis: how nonculture diagnostics will improve understanding of disease spectrum and transform patient care. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2013; **56**: 1284-92.
 39. Clancy CJ, Nguyen MH. Diagnosing Invasive Candidiasis. *J Clin Microbiol* 2018; **56**.
 40. Cockerill FR, 3rd, Wilson JW, Vetter EA et al. Optimal testing parameters for blood cultures. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2004; **38**: 1724-30.
 41. Collin J. Svårigheter i att detektera Candida albicans i blododling vid samtidig bakteriemi. . *Examensarbete BMA062 Sahlgrenska akademien* 2019.
 42. Connell TG, Rele M, Cowley D et al. How reliable is a negative blood culture result? Volume of blood submitted for culture in routine practice in a children's hospital. *Pediatrics* 2007; **119**: 891-6.
 43. Dawson S. Blood culture contaminants. *J Hosp Infect* 2014; **87**: 1-10.
 44. Dayer MJ, Jones S, Prendergast B et al. Incidence of infective endocarditis in England, 2000-13: a secular trend, interrupted time-series analysis. *Lancet* 2015; **385**: 1219-28.
 45. Diekema D, Arbefeville S, Boyken L et al. The changing epidemiology of healthcare-associated candidemia over three decades. *Diagnostic microbiology and infectious disease* 2012; **73**: 45-8.
 46. Dien Bard J, McElvania TeKippe E. Diagnosis of Bloodstream Infections in Children. *J Clin Microbiol* 2016; **54**: 1418-24.
 47. Dunne WM, Jr., Tillman J, Havens PL. Assessing the need for anaerobic medium for the recovery of clinically significant blood culture isolates in children. *The Pediatric infectious disease journal* 1994; **13**: 203-6.
 48. Durack DT, Lukes AS, Bright DK. New criteria for diagnosis of infective endocarditis: utilization of specific echocardiographic findings. Duke Endocarditis Service. *Am J Med* 1994; **96**: 200-9.
 49. Dwivedi S, Bhalla R, Hoover DR et al. Discarding the initial aliquot of blood does not reduce contamination rates in intravenous-catheter-drawn blood cultures. *J Clin Microbiol* 2009; **47**: 2950-1.

50. Edelmann A, Pietzcker T, Wellinghausen N. Comparison of direct disk diffusion and standard microtitre broth dilution susceptibility testing of blood culture isolates. *J Med Microbiol* 2007; **56**: 202-7.
51. Erb S, Frei R, Schregenberger K et al. Sonication for diagnosis of catheter-related infection is not better than traditional roll-plate culture: a prospective cohort study with 975 central venous catheters. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2014; **59**: 541-4.
52. Ericson EL, Klingspor L, Ullberg M et al. Clinical comparison of the Bactec Mycosis IC/F, BacT/Alert FA, and BacT/Alert FN blood culture vials for the detection of candidemia. *Diagnostic microbiology and infectious disease* 2012; **73**: 153-6.
53. ESCMID. European Manual of Clinical Microbiology. 1st edition. 2012.
54. Evans L, Rhodes A, Alhazzani W et al. Surviving sepsis campaign: international guidelines for management of sepsis and septic shock 2021. *Intensive care medicine* 2021; **47**: 1181-247.
55. Falagas ME, Ierodiakonou V, Alexiou VG. Clinical practice of obtaining blood cultures from patients with a central venous catheter in place: an international survey. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases* 2009; **15**: 683-6.
56. Fang H, Ohlsson AK, Ullberg M et al. Evaluation of species-specific PCR, Bruker MS, VITEK MS and the VITEK 2 system for the identification of clinical Enterococcus isolates. *European journal of clinical microbiology & infectious diseases : official publication of the European Society of Clinical Microbiology* 2012; **31**: 3073-7.
57. Fernández-Cruz A, Martín-Rabadán P, Suárez-Salas M et al. Is it feasible to diagnose catheter-related candidemia without catheter withdrawal? *Medical mycology* 2014; **52**: 491-7.
58. Folkhälsomyndigheten. Packa provet rätt. 2020; 22.
59. Gahlot R NC, Kumar V, Yadav G, Anupurba S. Catheter-related bloodstream infections. *Int J Crit Illn Inj Sci* 2014; **4**: 162-7.
60. Garcia RA, Spitzer ED, Beaudry J et al. Multidisciplinary team review of best practices for collection and handling of blood cultures to determine effective interventions for increasing the yield of true-positive bacteremias, reducing contamination, and eliminating false-positive central line-associated bloodstream infections. *Am J Infect Control* 2015; **43**: 1222-37.
61. Gonsalves WI, Cornish N, Moore M et al. Effects of volume and site of blood draw on blood culture results. *J Clin Microbiol* 2009; **47**: 3482-5.
62. Gowardman JR, Jeffries P, Lassig-Smith M et al. A comparative assessment of two conservative methods for the diagnosis of catheter-related infection in critically ill patients. *Intensive care medicine* 2013; **39**: 109-16.
63. Greub G, Lepidi H, Rovey C et al. Diagnosis of infectious endocarditis in patients undergoing valve surgery. *Am J Med* 2005; **118**: 230-8.
64. Gristina AG, Costerton JW. Bacterial adherence to biomaterials and tissue. The significance of its role in clinical sepsis. *J Bone Joint Surg Am* 1985; **67**: 264-73.
65. Guembe M, Martin-Rabadan P, Cruces R et al. Sonicating multi-lumen sliced catheter tips after the roll-plate technique improves the detection of catheter colonization in adults. *Journal of microbiological methods* 2016; **122**: 20-2.
66. Guembe M, Martin-Rabadan P, Cruces R et al. Corrigendum to "Sonicating multi-lumen sliced catheter tips after the roll-plate technique improves the detection of catheter colonization in adults" [J. Microbiol. Methods 122 (2016) 20-22]. *Journal of microbiological methods* 2016; **127**: 242.

67. Haimi-Cohen Y, Shafinoori S, Tucci V et al. Use of incubation time to detection in BACTEC 9240 to distinguish coagulase-negative staphylococcal contamination from infection in pediatric blood cultures. *The Pediatric infectious disease journal* 2003; **22**: 968-74.
68. Hall KK, Lyman JA. Updated review of blood culture contamination. *Clinical microbiology reviews* 2006; **19**: 788-802.
69. Hamed A, Mamoury GA, Akhlaghi F. A comparative study of blood culture sampling from umbilical catheter line versus peripheral site. *Acta medica Iranica* 2010; **48**: 231-3.
70. Hockey LJ, Fujita NK, Gibson TR et al. Detection of fungemia obscured by concomitant bacteremia: in vitro and in vivo studies. *J Clin Microbiol* 1982; **16**: 1080-5.
71. Hogevik H, Olaison L, Andersson R et al. Epidemiologic aspects of infective endocarditis in an urban population. A 5-year prospective study. *Medicine (Baltimore)* 1995; **74**: 324-39.
72. Hoiby N, Bjarnsholt T, Moser C et al. ESCMID guideline for the diagnosis and treatment of biofilm infections 2014. *Clinical microbiology and infection : the official publication of the European Society of Clinical Microbiology and Infectious Diseases* 2015; **21 Suppl 1**: S1-25.
73. Humphries RM, Kircher S, Ferrell A et al. The Continued Value of Disk Diffusion for Assessing Antimicrobial Susceptibility in Clinical Laboratories: Report from the Clinical and Laboratory Standards Institute Methods Development and Standardization Working Group. *J Clin Microbiol* 2018; **56**.
74. Isaacman DJ, Karasic RB, Reynolds EA et al. Effect of number of blood cultures and volume of blood on detection of bacteremia in children. *The Journal of pediatrics* 1996; **128**: 190-5.
75. Johnsson ATA, Wong AYW, Özenci V. The impact of delayed analysis of positive blood cultures on the performance of short-term culture followed by MALDI-TOF MS. *J Microbiol Methods* 2020; **177**: 106027.
76. Kaasch AJ, Rieg S, Hellmich M et al. Differential time to positivity is not predictive for central line-related Staphylococcus aureus bloodstream infection in routine clinical care. *J Infect* 2014; **68**: 58-61.
77. Kellogg JA, Manzella JP, Bankert DA. Frequency of low-level bacteremia in children from birth to fifteen years of age. *J Clin Microbiol* 2000; **38**: 2181-5.
78. Klingspor L, Lindback E, Ullberg M et al. Seven years of clinical experience with the Yeast Traffic Light PNA FISH: Assay performance and possible implications on antifungal therapy. *Mycoses* 2018; **61**: 179-85.
79. Klingspor L, Muhammed SA, Ozenci V. Comparison of the two blood culture systems, Bactec 9240 and BacT/Alert 3D, in the detection of Candida spp. and bacteria with polymicrobial sepsis. *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 2012; **31**: 2983-7.
80. Klingspor LH, J. Svampinfektioner (invasiva), diagnostik. Stockholm: Internetmedicin.se, [uppdaterad 2018-01-21]. <https://www.internetmedicin.se>
81. Klotz SA, Chasin BS, Powell B et al. Polymicrobial bloodstream infections involving Candida species: analysis of patients and review of the literature. *Diagnostic microbiology and infectious disease* 2007; **59**: 401-6.
82. Kumar A, Sharma RM, Jaideep CN et al. Diagnosis of central venous catheter-related bloodstream infection without catheter removal: A prospective observational study. *Med J Armed Forces India* 2014; **70**: 17-21.

83. Lalezari A, Cohen MJ, Svinik O et al. A simplified blood culture sampling protocol for reducing contamination and costs: a randomized controlled trial. *Clin Microbiol Infect* 2020; **26**: 470-4.
84. Lamy B, Dargère S, Arendrup MC et al. How to Optimize the Use of Blood Cultures for the Diagnosis of Bloodstream Infections? A State-of-the Art. *Front Microbiol* 2016; **7**: 697.
85. Lamy B, Roy P, Carret G et al. What is the relevance of obtaining multiple blood samples for culture? A comprehensive model to optimize the strategy for diagnosing bacteremia. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2002; **35**: 842-50.
86. Lancaster DP, Friedman DF, Chiotos K et al. Blood Volume Required for Detection of Low Levels and Ultralow Levels of Organisms Responsible for Neonatal Bacteremia by Use of Bactec Peds Plus/F, Plus Aerobic/F Medium, and the BD Bactec FX System: an In Vitro Study. *J Clin Microbiol* 2015; **53**: 3609-13.
87. Lee A, Mirrett S, Reller LB et al. Detection of bloodstream infections in adults: how many blood cultures are needed? *Journal of clinical microbiology* 2007; **45**: 3546-8.
88. Leyssene D, Gardes S, Vilquin P et al. Species-driven interpretation guidelines in case of a single-sampling strategy for blood culture. *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 2011; **30**: 1537-41.
89. Li J, Plorde JJ, Carlson LG. Effects of volume and periodicity on blood cultures. *Journal of clinical microbiology* 1994; **32**: 2829-31.
90. Li JS, Sexton DJ, Mick N et al. Proposed modifications to the Duke criteria for the diagnosis of infective endocarditis. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2000; **30**: 633-8.
91. Liesman RM, Pritt BS, Maleszewski JJ et al. Laboratory Diagnosis of Infective Endocarditis. *Journal of clinical microbiology* 2017; **55**: 2599-608.
92. Lin JN, Lai CH, Chen YH et al. Characteristics and outcomes of polymicrobial bloodstream infections in the emergency department: A matched case-control study. *Academic emergency medicine: official journal of the Society for Academic Emergency Medicine* 2010; **17**: 1072-9.
93. Liumbruno GM, Catalano L, Piccinini V et al. Reduction of the risk of bacterial contamination of blood components through diversion of the first part of the donation of blood and blood components. *Blood transfusion = Trasfusione del sangue* 2009; **7**: 86-93.
94. Lobmaier IV, Vege A, Gaustad P et al. Bacteriological investigation--significance of time lapse after death. *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 2009; **28**: 1191-8.
95. Läkemedelsverket. Neonatal sepsis – ny behandlingsrekommendation. Information från Läkemedelsverket. 2013;24(3):15–25.
96. Maki DG, Weise CE, Sarafin HW. A semiquantitative culture method for identifying intravenous-catheter-related infection. *N Engl J Med* 1977; **296**: 1305-9.
97. Melhus Å. Bakteriologisk diagnostik i neonatalperioden. Information från läkemedelsverket 2013;24(3):44-47.
98. Mellingshoff SC, Cornely OA, Jung N. Essentials in Candida bloodstream infection. *Infection* 2018; **46**: 897-9.
99. Menchinelli G, Liotti FM, Fiori B et al. In vitro Evaluation of BACT/ALERT® VIRTUO®, BACT/ALERT 3D®, and BACTEC™ FX Automated Blood Culture Systems for Detection of Microbial Pathogens Using Simulated Human Blood Samples. *Front Microbiol* 2019; **10**: 221.

100. Mermel LA, Maki DG. Detection of bacteremia in adults: consequences of culturing an inadequate volume of blood. *Ann Intern Med* 1993; **119**: 270-2.
101. Mitteregger D, Barousch W, Nehr M et al. Neutralization of antimicrobial substances in new BacT/Alert FA and FN Plus blood culture bottles. *J Clin Microbiol* 2013; **51**: 1534-40.
102. Monsen T, Lövgren E, Widerström M et al. In Vitro Effect of Ultrasound on Bacteria and Suggested Protocol for Sonication and Diagnosis of Prosthetic Infections. *Journal of Clinical Microbiology* 2009; **47**: 2496-501.
103. Morris AJ, Drinkovic D, Pottumarthy S et al. Gram stain, culture, and histopathological examination findings for heart valves removed because of infective endocarditis. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2003; **36**: 697-704.
104. Munoz P, Bouza E, Marin M et al. Heart valves should not be routinely cultured. *Journal of clinical microbiology* 2008; **46**: 2897-901.
105. Munson EL, Diekema DJ, Beekmann SE et al. Detection and treatment of bloodstream infection: laboratory reporting and antimicrobial management. *J Clin Microbiol* 2003; **41**: 495-7.
106. Noel GJ, Laufer DA, Edelson PJ. Anaerobic bacteremia in a neonatal intensive care unit: an eighteen-year experience. *The Pediatric infectious disease journal* 1988; **7**: 858-62.
107. Orihuela-Martín J, Rodríguez-Núñez O, Morata L et al. Performance of differential time to positivity as a routine diagnostic test for catheter-related bloodstream infections: a single-centre experience. *Clin Microbiol Infect* 2020; **26**: 383.e1-.e7.
108. Ortiz C, Lopez J, Garcia H et al. Clinical classification and prognosis of isolated right-sided infective endocarditis. *Medicine (Baltimore)* 2014; **93**: e137.
109. Ozenci V, Tegmark-Wisell K, Lundberg C et al. Rapid culture and identification: a practical method for early preliminary laboratory diagnosis of sepsis. *Clin Microbiol Infect* 2008; **14**: 177-80.
110. Pappas PG, Lionakis MS, Arendrup MC et al. Invasive candidiasis. *Nature reviews Disease primers* 2018; **4**: 18026.
111. Parienti JJ, du Cheyron D, Ramakers M et al. Alcoholic povidone-iodine to prevent central venous catheter colonization: A randomized unit-crossover study. *Critical care medicine* 2004; **32**: 708-13.
112. Patton RG, Schmitt T. Innovation for reducing blood culture contamination: initial specimen diversion technique. *J Clin Microbiol* 2010; **48**: 4501-3.
113. Pavlaki M, Poulakou G, Drimousis P et al. Polymicrobial bloodstream infections: Epidemiology and impact on mortality. *Journal of global antimicrobial resistance* 2013; **1**: 207-12.
114. Péan de Ponfilly G, Benmansour H, Manda V et al. Impact of 24/7 loading of blood culture bottles in a new automated incubator on the diagnosis of bloodstream infections. *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 2021.
115. Petti CA, Bhally HS, Weinstein MP et al. Utility of extended blood culture incubation for isolation of *Haemophilus*, *Actinobacillus*, *Cardiobacterium*, *Eikenella*, and *Kingella* organisms: a retrospective multicenter evaluation. *Journal of clinical microbiology* 2006; **44**: 257-9.
116. Quittnat Pelletier F, Joarder M, Poutanen SM et al. Evaluating Approaches for the Diagnosis of Hemodialysis Catheter-Related Bloodstream Infections. *Clinical journal of the American Society of Nephrology: CJASN* 2016; **11**: 847-54.

117. Raad I, Hanna H, Maki D. Intravascular catheter-related infections: advances in diagnosis, prevention, and management. *The Lancet Infectious diseases* 2007; **7**: 645-57.
118. Raad I, Hanna HA, Alakech B et al. Differential time to positivity: a useful method for diagnosing catheter-related bloodstream infections. *Ann Intern Med* 2004; **140**: 18-25.
119. Rand KH, Beal SG, Rivera K et al. Hourly Effect of Pretreatment With IV Antibiotics on Blood Culture Positivity Rate in Emergency Department Patients. *Open Forum Infect Dis* 2019; **6**: ofz179.
120. Razimi R. MP. ABC om Infektiös endokardit. *Läkartidningen* 2019; **116**: 1-5.
121. Riedel S. The value of postmortem microbiology cultures. *J Clin Microbiol* 2014; **52**: 1028-33.
122. Ronnberg C, Mildh M, Ullberg M et al. Transport time for blood culture bottles: underlying factors and its consequences. *Diagnostic microbiology and infectious disease* 2013; **76**: 286-90.
123. Roth A, Wiklund AE, Pålsson AS et al. Reducing blood culture contamination by a simple informational intervention. *J Clin Microbiol* 2010; **48**: 4552-8.
124. Rovey C, Greub G, Lepidi H et al. PCR detection of bacteria on cardiac valves of patients with treated bacterial endocarditis. *Journal of clinical microbiology* 2005; **43**: 163-7.
125. Rupp ME, Cavalieri RJ, Marolf C et al. Reduction in Blood Culture Contamination Through Use of Initial Specimen Diversion Device. *Clinical infectious diseases : an official publication of the Infectious Diseases Society of America* 2017; **65**: 201-5.
126. Ruus CS, S. Magnusson, T. Snygg-Martin, U. Studahl, M. Andersson, L-M. Tidig infektionskonsult gav effect vid Staphylococcus aureus-bakteriemi. *Läkartidningen* 2018; **3**: 88-91.
127. Sarkar S, Bhagat I, DeCristofaro JD et al. A study of the role of multiple site blood cultures in the evaluation of neonatal sepsis. *Journal of perinatology : official journal of the California Perinatal Association* 2006; **26**: 18-22.
128. Sautter RL, Bills AR, Lang DL et al. Effects of delayed-entry conditions on the recovery and detection of microorganisms from BacT/ALERT and BACTEC blood culture bottles. *J Clin Microbiol* 2006; **44**: 1245-9.
129. Savage TJ, Rao S, Joerger J et al. Predictive Value of Direct Disk Diffusion Testing from Positive Blood Cultures in a Children's Hospital and Its Utility in Antimicrobial Stewardship. *J Clin Microbiol* 2021; **59**.
130. Scheer CS, Fuchs C, Grundling M et al. Impact of antibiotic administration on blood culture positivity at the beginning of sepsis: a prospective clinical cohort study. *Clin Microbiol Infect* 2019; **25**: 326-31.
131. Schelonka RL, Chai MK, Yoder BA et al. Volume of blood required to detect common neonatal pathogens. *The Journal of pediatrics* 1996; **129**: 275-8.
132. Sherertz RJ, Raad, II, Belani A et al. Three-year experience with sonicated vascular catheter cultures in a clinical microbiology laboratory. *Journal of clinical microbiology* 1990; **28**: 76-82.
133. Silver H SA. A practical and efficacious method for obtaining significant postmortem blood cultures. *American journal of clinical pathology* 1969; **52**: 433-7.
134. Sjölin J HL. Svampinfektioner (invasiva), behandling. Göteborg. [uppdaterad 2021-10-12]. <https://www.internetmedicin.se>.
135. Stempel JM, Farmakiotis D, Tarrand JJ et al. Time-to-reporting of blood culture positivity and central venous catheter-associated *Candida glabrata* fungemia in cancer patients. *Diagnostic microbiology and infectious disease* 2016; **85**: 391-3.

136. Syed S, Liss DT, Costas CO et al. Diversion Principle Reduces Skin Flora Contamination Rates in a Community Hospital. *Archives of pathology & laboratory medicine* 2020; **144**: 215-20.
137. Taniguchi T, Tsuha S, Shiiki S et al. High positivity of blood cultures obtained within two hours after shaking chills. *International journal of infectious diseases: IJID : official publication of the International Society for Infectious Diseases* 2018; **76**: 23-8.
138. Towns ML, Jarvis WR, Hsueh PR. Guidelines on blood cultures. *J Microbiol Immunol Infect* 2010; **43**: 347-9.
139. Venturelli C, Righi E, Borsari L et al. Impact of Pre-Analytical Time on the Recovery of Pathogens from Blood Cultures: Results from a Large Retrospective Survey. *PloS one* 2017; **12**: e0169466.
140. Vondracek M, Sartipy U, Aufwerber E et al. 16S rDNA sequencing of valve tissue improves microbiological diagnosis in surgically treated patients with infective endocarditis. *The Journal of infection* 2011; **62**: 472-8.
141. Werner M, Andersson R, Olaison L et al. A 10-year survey of blood culture negative endocarditis in Sweden: aminoglycoside therapy is important for survival. *Scand J Infect Dis* 2008; **40**: 279-85.
142. Yu D, Larsson A, Parke Å et al. Single-Sampling Strategy vs. Multi-Sampling Strategy for Blood Cultures in Sepsis: A Prospective Non-inferiority Study. *Front Microbiol* 2020; **11**: 1639.
143. Zaidi AK, Knaut AL, Mirrett S et al. Value of routine anaerobic blood cultures for pediatric patients. *The Journal of pediatrics* 1995; **127**: 263-8.
144. Zimmerman FS, Karamah H, Ben-Chetrit E et al. Modification of blood test draw order to reduce blood culture contamination: a randomized clinical trial. *Clinical infectious diseases: an official publication of the Infectious Diseases Society of America* 2019.